

DNA Isothermal Amplification Kit (Basic Version)

Product Name

Common Name: DNA Isothermal Amplification Kit (Basic Version)

Specification

- **Item No.:** WLB8201KIT
- **Specification:** 48 rxns/kit

Description

The DNA Isothermal Amplification Kit (Basic Version) is developed based on a room-temperature isothermal nucleic acid rapid amplification technology. It is suitable for laboratory-grade DNA amplification and other detection applications requiring DNA amplification.

Principal

This kit is based on rapid nucleic acid amplification technology that operates at room temperature: under isothermal conditions at room temperature, the recombinase and primers form a Rec/ssDNA complex (protein/single-stranded nucleic acid). With the assistance of accessory proteins and single-stranded binding proteins (SSB), this complex invades the double-stranded DNA template. A D-loop region forms at the invasion site, initiating scanning of the DNA duplex. Once the target region complementary to the primer is identified, the Rec/ssDNA complex disassembles, allowing the polymerase to bind to the 3' end of the primer and initiate extension.

Features and Advantages

- This kit offers high sensitivity, strong specificity, and a rapid reaction time (only takes 30 minutes). The reagents are provided in a lyophilized form for easy handling and storage.
- Operates with basic equipment such as metal bath or water bath, without the need for expensive

equipment like PCR machines.

Primer Design

Primer with a length of 30-35 bp is recommended, since too short primer can reduce amplification speed and detection sensitivity. Formation of secondary structures should be avoided during design of primer to ensure the amplification efficient. The length of amplicon is suggested to be 150-300bp, typically not exceeding 500bp.

Storage

- **Shipping Temperature:** Stable temperature of ≤ 20 °C.
- **Storage Conditions:** Store at ≤ -20 °C (± 5 °C) in a constant temperature environment, protected from light, and avoid heavy pressure and repeated freeze-thaw cycles.
- **Shelf Life:** 14 months.
- **Production Date:** See the outer packaging.

Kit Component

Component	Volume
A buffer	1.6 mL×1 vial
B buffer	150 μ L×1 vial
Positive control template	30 μ L×1 vial
Positive control primer mix	60 μ L×1 vial
Total	48 rxns
User manual	1 copy

Procedure

Thaw the necessary reagents of the kit at room temperature 30 minutes in advance. Vortex to mix thoroughly.

1. Add 29.4 μ L of A buffer to each lyophilized reaction tube.

Note: Ensure the A buffer thawed and mixed thoroughly, as incomplete mixing may affect the

experimental results.

2. Add 2 μL of forward primer and 2 μL of reverse primer with a concentration of 10 μM to each reaction tube.

Note: For multiple reactions, combine steps 1 and 2, then aliquot into the respective tubes.

3. Add 2-14.1 μL of nucleic acid template to the reaction tubes, the volume can be adjusted as needed. Correspondingly adjust the volume of ddH₂O to ensure the total volume of the template and ddH₂O is 14.1 μL .

4. Add 2.5 μL of B buffer to each reaction tube and mix thoroughly.

Note: Invert the tube 8-10 times to ensure thorough mixing. For multiple reactions, it is recommended to add the B buffer to the inner side of the reaction tube cap, then close the tube and mix by inverting.

5. After mixing, briefly spin down the reaction mixture to the bottom of the tube, or perform a quick centrifugation. Immediately and gently place the reaction tubes in an incubator set at 37-39 °C for 30 minutes.

6. After the reaction, add an equal volume of Tris-saturated phenol/chloroform/isoamyl alcohol in a ratio of 25:24:1 to the reaction mixture. Mix thoroughly in a one-to-one ratio, then centrifuge at 12,000 rpm for 5 minutes. Take 5 μL of the supernatant for agarose gel electrophoresis with a recommended agarose gel concentration of 1.5-2%.

Note:

- High-temperature denaturation may be ineffective at removing proteins, potentially impacting the analysis of results.
- Some commercial kits for purifying amplification products may not be compatible, which could result in false-negative outcomes.

PCR Configuration


Reagent	Volume (μL)
A buffer	29.4
forward primer (10 μM)*	2

reverse primer (10 μM)*	2
ddH ₂ O and DNA template	14.1
B buffer	2.5
Total Volume	50

*** Preparation of Positive Control Reaction Mix**

Add 2 μL of the positive control template and 4μL of the positive control primer mix (including both forward and reverse primers). Prepare the remaining components according to the protocol above.

Positive Control Electrophoresis Results



M: Marker P: Positive Control N: Negative Control

Precautions

- Given the high sensitivity of the kit, it is essential to avoid nucleic acid contamination during the reaction and to include a blank control.
- Only remove the necessary number of MIRA reaction units for your experiment, and store any unused units under the recommended conditions.